

# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination of Isolator Systems*

Volker Sigwarth  
Skan AG  
Switzerland



## *Process Development*

# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination Process Development*

- Overview of current Regulations and Standards
- *Biological Indicators*
- *Investigation on Germ reducing Effects*
- *Method for Cycle Development*
- Discussion of Risks, Possibilities, and Experience

# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination Process Development*

- Research Article                      PDA Journal

“Effect of Carrier Materials on the Resistance of Spores of *Bacillus Stearothermophilus* to gaseous Hydrogen Peroxide”

Volker Sigwarth, Skan AG

Alexandra Stärk, Novartis Pharma AG

PDA Journal, Vol. 57, No.1 January / February 2003

# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination Process Development*

- *FDA*; GMP Guidance for Industry; Draft August 2003

*“Sterile Drug Products Produced by Aseptic Processing”*

- *Decontamination Efficacy*

“Process development and validation studies should include a thorough determination of cycle capability. The characteristics of these agents generally preclude the reliable use of statistical methods (e.g. fractional negative) to determine the process lethality (Ref. 13).”

# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination Process Development*

- **Research Article**                      **PDA Journal**

“Effect of Carrier Materials on the Resistance of Spores of *Bacillus Stearothermophilus* to gaseous Hydrogen Peroxide”

Volker Sigwarth, Skan AG

Alexandra Stärk, Novartis Pharma AG

PDA Journal, Vol. 57, No.1 January / February 2003

- **Cited as Reference (13);**              **FDA GMP Guidance for Industry**

# *Overview of current Regulation and Standards*

## *Process Development*


- EN / ISO 14161 “Guidance for Biological Indicators”
- USP <1208> „Sterility Testing - Validation of Isolator Systems”  
<55> “ Biological Indicators - Resistance Performance Test”  
<1035> “Biological Indicators for Sterilization”
- FDA “Sterile Drug Produced by Aseptic Processing”  
Appendix 1: “Aseptic Processing Isolators”
- PIC/S Isolators Used For Aseptic Processing and Sterility Testing

# *Overview of current Regulation and Standards Requirements for Validation*

- Definition of Process- *Performance* and *Boundaries*
- Proof of Process- *Performance* within *Boundaries*
- Control of Process- *Parameters* and *Influences*
  
- *Process Comprehension for individual Application*

# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination Process Development*

- Process Development of an alternative Sterilization Method





# *Biological Indicators for alternative Sterilization*

- H<sub>2</sub>O<sub>2</sub> Decontamination as alternative Sterilization Method
- Model of microbial Reduction                      Biological Indicators
- D-value Determination                              Biological Indicators
- Composition of Biological Indicators  
    Test Organism; Initial Population; Carrier Material; Primary Packaging
- Samples commercial available Biological Indicators



## *Alternative Sterilization Methods*

- Chemical Processes
  - Ethylene Oxide
  - Per Acetic Acid
  - Hydrogen Peroxide, H<sub>2</sub>O<sub>2</sub>
  - Ozone
  - Chlorine Dioxide
- Physical Processes
  - Electron Beam
  - Micro Waves
  - UV Beam
- *Use only possible if Standard Methods are not applicable*

## *Alternative Sterilization Methods*


- One Step      *lower*      than Standard Sterilization Methods
- Terminology      Sanitization  
                                 Decontamination  
                                 Inactivation  
                                 Disinfection
- Process Result      „Spore Log Reduction“      *SLR*
- *Based on the Process Expectations      individually defined*

# *Alternative Sterilization Methods*

## *H<sub>2</sub>O<sub>2</sub> Decontamination of Isolator Systems*

- Surface Decontamination of the Isolator Chamber
- Vaporizing of aqueous H<sub>2</sub>O<sub>2</sub> Solution
- Sporicidal Inactivation Process
- widely used in pharmaceutical Industry
  
- *Process Parameters are individually applied*

# Isolator System



Material: stainless steel, glass

Volume: 1,4 m<sup>3</sup> (40 ft<sup>3</sup>)

# *Process Control*


## *Process*

- Temperature [°C]
- Humidity [% rH]
- Air Velocity [m/s]
- Pressure [Pa]
- Mass / Balance [g]

## *H<sub>2</sub>O<sub>2</sub> Gas-Concentration*


- Electro-chemical Sensor
- UV-Spectrometer
- IMS-Spectrometer
- NIR-Spectrometer
- Wet-chemistry Method

# *Decontamination Cycle*






# *Decontamination Cycle*




Phase 1:            Pre-conditioning  
                         to establish the initial conditions in the chamber

# *Decontamination Cycle*




Phase 2:            Conditioning  
                          to establish the decontamination effect

# *Decontamination Cycle*



Phase 3:            **Decontamination**  
to maintain stable decontamination effect  
to ensure the total bacterial reduction over time

# *Decontamination Cycle*




Phase 4:            **Aeration**  
to reach minimal residual H<sub>2</sub>O<sub>2</sub> concentration

## *Alternative Sterilization Methods*

- No useful Correlation of:

*Physical Parameters*




*Inactivation Effect*

- Relevant Process Parameters and Values                      not known
- Range and Boundaries of Process Values                      not known
- Design, Qualification, Routine Use                              generates Problems
- Unexpected high Effort for the Qualification of a single Application

## *Alternative Sterilization Methods*

- No useful Correlation of:

*Physical Parameters*



*Inactivation Effect*

- The use of alternative Inactivation Methods requires often:
- *Description of the Inactivation Effect directly*
- *Measure of the Inactivation Effect*
- *by the Use of Biological Indicators*

# *Microbiological System*

## *Description of Biological Indicator*

- defined Test Organism
- defined initial Population
- Carrier Material
- Primary Packaging
- defined Resistance to a specified Inactivation Method

# *Microbiological System*

## *Resistance Description of Biological Indicator*


- Initial Population of the Test Organisms [CFU/Carrier]
- D-value [min]
- Survival - Kill Window [min]
  - Survival time [min] = D-value x ([log] Population – 2)
  - Kill time [min] = D-value x ([log] Population + 4)




# Microbiological System

## Model of Microbial Reduction

- Initial Population [log-scale]
- Inactivation Time [min]
- Survival Curve
- D-value [min]
- Survival - Kill Window [min]




# Microbiological System



# Microbiological System

## Model of Microbial Reduction

- Initial Population [log-scale]
- Inactivation Time [min]
- Survival Curve
- D-value [min]
- Survival - Kill Window [min]



# *Microbiological System*

## D-value of Biological Indicators

- *Measure* of the Inactivation Effect
- *Quantification* achieved Inactivation Effect

## *Important*

- D-value Slope of the Survival Curve
- corresponding Model Behavior survival, fractional, kill

## *Methods of D-value Determination*

- Numeration of the residual Population after Inactivation

Survivor Curve Method

*SCM*

- Fractional Negative Methods

Stumbo Murphy Cochran Method

*SMCM*

Limited Spearman Karber Method

*LSKM*

*Suitable Method*


Selection of Biological Indicators

Description of the Inactivation Effect

# *D-value Determination*

## *SC Method*

- defined numbers of BI
- defined Exposure Times
- Numeration of residual Population
- Plot of Survivor Curve
  
- < 50% of initial Population
- > 50 Counts / Carrier




## *D-value Determination*

### *SC Method*

- Survival Window
- Beginning of fractional Field
- high Lab Effort
- good if *Total Kill* is not required
- Disinfection Testing
- “*Last Chance Method*”

No Information about:


- *Kill Window*
- *Model Behaviour*



## *D-value Determination*

### *SMCM*

- high Number of BIs
- one Exposure Time; fractional Field
- Relationship of survival / kill
- initial Population
- Calculate D-value






## *D-value Determination*

### *SMC Method*

- initial Population
- Survival Window
- Beginning of fractional Field
- high statistical Accuracy
- low Lab Effort
- D-value has to be known
- good if *Total Kill* is not required

No Information about:

- *Kill Window*
- *Model Behaviour*



# Stumbo Murphy Cochran Method

# SMCM

Initial Population  $N_0: 1.0 \times 10^6$

exposure	01										
exposure time [min]	13										
result	-	-	-	-	+	+	+	+	+	+	neg
	-	-	-	-	+	+	+	+	+	+	-
	-	-	-	+	+	+	+	+	+	+	pos
	-	-	-	+	+	+	+	+	+	+	+
	-	-	-	+	+	+	+	+	+	+	


+ growth  
- no growth

D-value = 2.03 [min]  
95% CI D-value = 2.03 ± 0.05 [min]

# D-value Determination

## LSKM


- several Groups of Bioindicators
- different Exposure Times
- observed Model Behavior
- initial Population
- D-value Calculation



# D-value Determination

## LSKM

- Kill Window is covered
- Model Behavior is completely shown



# Limited Spearman Karber Method LSKM

Initial Population

$N_0: 1.0 \times 10^6$

exposure	01	02	03	04	05	06	07	08	09	10	
exposure time [min]	6.0	8.5	11.0	13.5	16.0	18.5	21.0	23.5	26.0	28.5	
result	1	+	+	+	+	+	-	-	-	-	neg
	2	+	+	+	+	+	-	-	-	-	-
	3	+	+	+	+	-	-	-	-	-	-
	4	+	+	+	+	-	-	-	-	-	pos
	5	+	+	+	-	-	-	-	-	-	+
	6	+	+	+	-	-	-	-	-	-	+
	7	+	+	+	-	-	-	-	-	-	
	8	+	+	+	-	-	-	-	-	-	
	9	+	+	-	-	-	-	-	-	-	
	10	+	-	-	-	-	-	-	-	-	

D-value = 2.06 [min]

95% CI D-value =  $2.06 \pm 0.22$  [min]

+ growth  
- no growth

# *Methods of D-value Determination*

## *LSKM*

- Determination the *Resistance* of Biological Indicator
- Evaluation of the *complete Model Behaviour*
  
- Certification of commercial Biological Indicators
- High Effort and Costs

# *Description of the Inactivation Effect*

Selection of Biological Indicators

Description of the Inactivation Effect

- good Estimation of the D- value
- good Understanding of the Model Behavior
- useful and pragmatic Tool
- good Relationship between Information and Costs

# Limited Spearman Karber Method LSKM

Initial Population

$N_0: 1.0 \times 10^6$

exposure	01	02	03	04	05	06	07	08	09	10	
exposure time [min]	6.0	8.5	11.0	13.5	16.0	18.5	21.0	23.5	26.0	28.5	
result	1	+	+	+	+	+	-	-	-	-	neg
	2	+	+	+	+	+	-	-	-	-	-
	3	+	+	+	+	-	-	-	-	-	-
	4	+	+	+	+	-	-	-	-	-	pos
	5	+	+	+	-	-	-	-	-	-	+
	6	+	+	+	-	-	-	-	-	-	+
	7	+	+	+	-	-	-	-	-	-	
	8	+	+	+	-	-	-	-	-	-	
	9	+	+	-	-	-	-	-	-	-	
	10	+	-	-	-	-	-	-	-	-	

D-value = 2.06 [min]

95% CI D-value =  $2.06 \pm 0.22$  [min]

+ growth  
- no growth



# Minimized LSKM, Reactive Pattern

Initial Population  $N_0: 1.0 \times 10^6$

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time[min]	6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
Result 1	+	+	+	+	-	-	-	-	-	-	neg
2	+	+	+	-	-	-	-	-	-	-	-
3	+	+	-	-	-	-	-	-	-	-	

+ growth  
- no growth

estimated D-value = 2.0 [min]

- statistical Accuracy lower
- lost of Information in the fractional Field
- good Estimation of the D-value

# *Description of the Inactivation Effect*

Selection of Biological Indicators

Description of the Inactivation Effect

- based on the Application of the minimized LSKM
- achieved Inactivation Effect can be *quantified*
- *Model Behavior* of Biological Indicator can be evaluated
- Application is useful and quite easy
- Good Relationship between Information and Costs

## *Selection of Biological Indicators*

- Test Organism
- Initial Population
- Carrier Material
- Primary Packaging

*Composition of the Biological Indicator has to reflect the Process Expectations*

## *Process Expectations*

### *H<sub>2</sub>O<sub>2</sub> Decontamination of Isolator Systems*

- Reduction of the microbial Contamination of the Isolator Chamber
- Reduction of the microbial Contamination on *Surfaces*

#### *No Penetration of the Inactivation Effect*

*FDA:* “Decontamination can be accomplished using a number of vaporizing agents, although these agents possess limited capability to penetrate obstructed or covered surfaces”

## *Process Expectations*

### *H<sub>2</sub>O<sub>2</sub> Decontamination of Isolator Systems*

- Reduction of the microbial Contamination of the Isolator Chamber
- Reduction of the microbial Contamination on *Surfaces*

#### *No Penetration of the Inactivation Effect*

- Total Kill of a 6 log Population *10 log Reduction*

# *Selection of Bioindicators*      *Reference Isolator*

- defined and controlled    *Reference Isolator*
- described and proven    *H<sub>2</sub>O<sub>2</sub> Decontamination Cycle*  
achieved Decontamination Effect  
Stability of Decontamination Effect

*Comparability of all following data*

## *Selection of Bioindicators*

## *Test Organism*

### *Pharmacopoeia USP < 1035 >*

- Test Organism                    „*highly resistance*”
- *vegetative Microbes*
- *Bacteria Spores*

*PIC/S:* “An understanding of the relationship between the resistance of the bioburden and that of the BI should be developed from trials and/or the literature”

## *Selection of Bioindicators*

## *Test Organism*

### *Vegetative Microbes*

Test Organism	D-value [min]	Carrier Material
Staphylococcus aureus	< 0.3	Glass
Pseudomonas aeruginosa	< 0.3	Glass
Micrococcus luteus	0.3 – 0.6	Glass
Acinetobacter lwoffii	0.6	Glass
Penicillium species	0.3 – 0.6	Glass
Aspergillus species	0.3	Glass

Method Deviation  $\pm$  10%



## *Selection of Bioindicators*

## *Test Organism*

### *Bacteria Spores*

Test Organism	D-value [min]	Carrier Material
Bacillus sphaericus	1.2	Glass
Bacillus subtilis var. niger	1.3	Glass
Bacillus stearothermophilus ATCC 12980	1.3	Glass
Bacillus stearothermophilus ATCC 7953	1.4	Glass

Method Deviation  $\pm$  10%

## *Selection of Bioindicators*

## *Test Organism*

- vegetative Microbes      *significant lower Resistance than B. Spores*
- Bacteria Spores      *comparable Resistance*

## *Selection of Bioindicators*

## *Test Organism*

### *Bacillus stearothermophilus*

- generally applied and accepted
- highly stable
- Incubation Temperature 55-60 °C
- Selective against Cross Contamination                      aseptic Handling

- *ATCC Strain 12980*

*generally used in the USA*

- *ATCC Strain 7953*

*generally used in Europe*

# Biological Indicator for gaseous H<sub>2</sub>O<sub>2</sub>

Testorganism: B. stearothermophilus ATCC 12980, min 1.0x10<sup>6</sup>  
 Carrier / Package: Stainless steel / Tyvek  
 Spezified D-value: 0.9 to 1.8 [min]

+ growth  
 - no growth

## Lot 01

exposure		01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]		6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
result	1	+	+	+	+	-	-	-	-	-	-	neg
	2	+	+	+	-	-	-	-	-	-	-	-
	3	+	+	-	-	-	-	-	-	-	-	

## Lot 02

exposure		01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]		6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
result	1	+	+	-	-	-	-	-	-	-	-	neg
	2	+	-	-	-	-	-	-	-	-	-	-
	3	+	-	-	-	-	-	-	-	-	-	

# *Biological Indicator for gaseous H<sub>2</sub>O<sub>2</sub>*

Testorganism: B. stearothermophilus ATCC 7953, min 1.0x10<sup>6</sup>  
 Carrier / Package: Stainless steel / Tyvek  
 Spezified D-value: 1.0 to 1.7 [min]

+ growth  
 - no growth

## Lot 01

exposure		01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]		6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
result	1	+	+	+	+	-	+	-	-	-	-	neg
	2	+	+	+	+	-	-	-	-	-	-	-
	3	+	+	+	-	-	-	-	-	-	-	

## Lot 02

exposure		01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]		6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
result	1	+	-	-	-	-	-	-	-	-	-	neg
	2	+	-	-	-	-	-	-	-	-	-	-
	3	+	-	-	-	-	-	-	-	-	-	

## *Selection of Bioindicators*

## *Test Organism*

### *Bacillus stearothermophilus*

commercially available Bioindicators with identical Composition

- bigger Differences of Resistance between *single Lots*
- than between *different Strains*

## *Selection of Bioindicators*

## *Population, $N_0$*


- Initial Population                      „*total kill*“ of a 6 log Population ?
- D-value                      *independent*                      of the initial Population
- D-value                      *only depends*                      Resistance of the individual Microbe

*Assumptions based on the ideal Model      empiric*

## *Initial Population, $N_0$*

### *Model of microbial Reduction*

- Survival Curve is a *straight Line*
- BI Resistance *independent* of  $N_0$
- *individual* Resistance of Microbe defines D-value of BI`s
- *Slope* of Survival Curve identically observed





## *Selection of Bioindicators*

*Population,  $N_0$*

### *Resistance in dependency of Initial Population*


Test Organism	Population, $N_0$ [log Steps]	D-value [min]
Micrococcus luteus	$\geq 1.0 \times 10^3$	0.3
Micrococcus luteus	$\geq 1.0 \times 10^4$	0.3
Micrococcus luteus	$\geq 1.0 \times 10^5$	0.6
Micrococcus luteus	$\geq 1.0 \times 10^6$	> 3.6

Method Deviation  $\pm 10\%$

## *Initial Population, $N_0$*

## *Model of microbial Reduction*

- Resistance of BI seems to *depend* on  $N_0$
- Survival Curve seems to be *not straight*
- Model becomes *questioned*



# Population $N_0$      Model Behavior

Test Organism:      B. stearothermophilus  
 carrier Material:      Glass

+ growth  
 - no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time [min]	3.0	6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	+
Result	1	+	+	+	-	-	-	-	-	-	neg
	2	+	+	-	-	-	-	-	-	-	-
	3	+	-	-	-	-	-	-	-	-	

- Initial Population:       $1 \times 10^5$
- estimated D-value:      1.3 [min]
- Model behavior:      OK

# Population $N_0$      Model Behavior

Test Organism:      B. stearothermophilus  
 carrier Material:      Glass

+ growth  
 - no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure time [min]	3.0	6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	+
Result 1	+	+	+	-	-	-	+	-	+	-	neg
2	+	+	+	-	-	-	+	-	-	-	-
3	+	+	-	-	-	-	-	-	-	-	

- Initial Population:       $1 \times 10^6$
- estimated D-value:      1.9 [min]
- Model Behavior:      not OK

## *Selection of Bioindicators*

## *Population, $N_0$*

By Increasing the Initial Population of BI's it can be observed

- *Increase in Resistance* higher D-value
- *Artifacts in the Model Behavior* late positive

H<sub>2</sub>O<sub>2</sub> Decontamination Process seems to be *not confidence*

## *Selection of Bioindicators*


## *Population, $N_0$*

### Dependency of Initial *Population* and *Resistance*

- vegetative Microbes are relatively *big* forming Agglomerates
- Preparation significantly *difficult* forming Agglomerates
- Suspensions often highly *dirty* forming a Coating


## *Selection of Bioindicators*

*Population,  $N_0$*



## *Selection of Bioindicators*


*Population,  $N_0$*





# *Selection of Bioindicators*

*Population,  $N_0$*



## *Selection of Bioindicators*

## *Population, $N_0$*

### Dependency of Initial *Population* and *Resistance*

- vegetative Microbes are relatively *big* forming Agglomerates
- Preparation significantly *difficult* forming Agglomerates
- Suspensions often highly *dirty* forming a Coating
- Agglomerates and Coating are *not penetrated* by  $H_2O_2$
- *inconsistent* Results of Bioindicators

The *Penetration* of the Decontamination Effect of the  $H_2O_2$  Process is

- *not ensured*
- *but also not expected*

# *Selection of Bioindicators*

# *Population, $N_0$*

## *Dependency of Initial Population and Resistance*

Test Organism	Population, $N_0$ [log Steps]	D-value [min]
Bacillus stearothermophilus	$\geq 1.0 \times 10^2$	1.6
Bacillus stearothermophilus	$\geq 1.0 \times 10^4$	1.3
Bacillus stearothermophilus	$\geq 1.0 \times 10^6$	1.4

Method Deviation  $\pm 10\%$


*Selection of Bioindicators*

*Population,  $N_0$*




# *Selection of Bioindicators*

*Population,  $N_0$*



# *Selection of Bioindicators*

*Population,  $N_0$*



## *Selection of Bioindicators*

## *Carrier Material*

### *PIC/S Recommendation on Isolator Technology*

“The carrier type e.g. plastic, paper, metal or other, of the biological indicator organism should be relevant to the materials being gassed or shown to be irrelevant”

### *FDA Comments on Isolator Technology*

“Rationale and justification of the use of stainless steel coupons as the challenge carrier for biological indicators”

## *Selection of Bioindicators*

## *Carrier Material*

### *Isolator Systems*

- Stainless- Steel      Isolator Chamber, Filling line, dif. Equipment
- Glass      Isolator Windows, div. Product Enclosures
- dif. Plastic Materials      Gloves, Sleeves, div. Product Enclosures
- dif. Metals      e.g. Aluminum

### *Biological Indicator*

- Should *reflect* the Process
- Should behave according the *Model of microbial Reduction*



# *Effect of Carrier Materials on the Resistance of B.stearothermophilus to gaseous H<sub>2</sub>O<sub>2</sub>*

- Research Article PDA Journal

“Effect of Carrier Materials on the Resistance of Spores of Bacillus Stearothermophilus to gaseous Hydrogen Peroxide”

Volker Sigwarth, Skan AG

Alexandra Stärk, Novartis Pharma AG

PDA Journal, Vol. 57, No.1 January / February 2003

# *Selection of Bioindicators*

# *Carrier Material*

Materials	Use
CrNi steel, different qualities: 1.4301 (304) unpolished 1.4301 (304) polished (Ra < 0.8 µm) 1.4435 (316L) unpolished 1.4435 (316L) polished (Ra < 0.8 µm)	Main chamber material Parts of filling line Steritest pump
Glass	Window and door material Media bottles, product units
Polycarbonate, PC	Window material
Hypalon	Glove material
Polyvinylchloride, PVC, soft	Material of glove gauntlets
Polyvinylchloride, PVC, hard	Glove ports
Polyvinylchloride, PVC	Package of steritest units
Polyoxymethylene, POM	Conveyer system of filling line, shelves of steritest isolator
Polypropylene, PP	RTP-ports
Polyethylene, PE-UHMW	Conveyer system of filling line, RTP-ports
Polytetrafluorethylene PTFE / (Teflon)	Tubings, parts of filling line
Aluminium anodized, different qualities	Material sample commercially available Material sample filling line Material sample air sampler
Butyl caoutchouc	Stoppers of media bottles and product units
Laminated foil 1; from inside to outside: Polyethylene, Aluminium, Polyester	Package of media plates for cleanroom monitoring
Laminated foil 2; from inside to outside: Polyethylene, Polypropylene	Package of media plates for cleanroom monitoring
Tyvek	Package of Steritest units
HEPA-filter pad	HEPA-filter


# *Carrier Material; Glass*

## Roughness

- P-P: 137 nm
- RMS: 3 nm
- $R_a$ : 1 nm

## Wettability

- high



# Carrier Material, Model Behavior

Test Organism: B. stearothermophilus  
Initial Population: > 1.0 x 10<sup>6</sup>  
Carrier Material: Glass

+ growth  
- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time [min]	2.0	4.0	6.0	8.0	10.0	12.0	14.0	16.0	18.0	20.0	+
Result 1	+	+	+	+	-	-	-	-	-	-	neg
2	+	+	+	-	-	-	-	-	-	-	-
3	+	+	+	-	-	-	-	-	-	-	

- estimated D-value: 1.2 [min]
- Model Behavior: OK


## *Carrier Material; stainless steel 1.4435, not polished*

### Roughness

- P-P: 817 nm
- RMS: 82 nm
- $R_a$ : 57 nm

### Wettability

- high



# Carrier Material, Model Behavior

Test Organism: B. stearothermophilus  
Initial Population: > 1.0 x 10<sup>6</sup>  
Carrier Material: Stainless Steel, not polished

+ growth  
- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time [min]	2.0	4.0	6.0	8.0	10.0	12.0	14.0	16.0	18.0	20.0	+
Result 1	+	+	+	+	+	-	-	-	-	-	neg
2	+	+	+	-	-	-	-	-	-	-	-
3	+	+	-	-	-	-	-	-	-	-	

- estimated D-value: 1.3 [min]
- Model Behavior: OK


## *Carrier Material; stainless steel 1.4435, polished*

### Roughness

- P-P: 1064 nm
- RMS: 165 nm
- $R_a$ : 127 nm

### Wettability

- high



# Carrier Material, Model Behavior

Test Organism: B. stearothermophilus  
Initial Population: > 1.0 x 10<sup>6</sup>  
Carrier Material: Stainless Steel, polished

+ growth  
- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Expousure [min]	2.5	5.0	7.5	10.0	12.5	15.0	17.5	20.0	22.5	25.0	+
Result 1	+	+	+	-	-	-	-	-	-	-	neg
2	+	+	-	-	-	-	-	-	-	-	-
3	+	-	-	-	-	-	-	-	-	-	

- estimated D-value: 0.9 [min]
- Model Behavior: OK




## *Carrier Material; PTFE, Teflon*

### Roughness

- P-P: 1266 nm
- RMS: 180 nm
- $R_a$ : 138 nm

### Wettability

- low



# Carrier Material, Model Behavior

Test Organism: B. stearothermophilus  
Initial Population: > 1.0 x 10<sup>6</sup>  
Carrier Material: PTFE

+ growth  
- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Expousure [min]	2.5	5.0	7.5	10.0	12.5	15.0	17.5	20.0	22.5	25.0	+
Result 1	+	+	+	+	+	-	-	-	-	-	neg
2	+	+	+	+	-	-	-	-	-	-	-
3	+	+	+	-	-	-	-	-	-	-	

- estimated D-value: 1.6 [min]
- Model Behavior: OK


# *Carrier Material, Hypalon*

## Roughness

- P-P: 2114 nm
- RMS: 350 nm
- $R_a$ : 286 nm

## Wettability

- low



# Carrier Material, Model Behavior

Test Organism: B. stearothermophilus  
Initial Population: > 1.0 x 10<sup>6</sup>  
Carrier Material: Hypalon

+ growth  
- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Expousure [min]	10.0	20.0	30.0	40.0	50.0	60.0	70.0	80.0	90.0	100.0	+
Result 1	+	+	-	-	-	-	-	-	-	-	neg
2	+	-	-	-	-	-	-	-	-	-	-
3	+	-	-	-	-	-	-	-	-	-	

- estimated D-value: 3.0 [min]
- Model Behavior: OK

## *Carrier Material,*


## *Aluminium; commercial Sample*

### Roughness

- P-P: 1322 nm
- RMS: 188 nm
- $R_a$ : 145 nm

### Wettability

- high



# Carrier Material, Model Behavior

Test Organism: B. stearothermophilus  
Initial Population: > 1.0 x 10<sup>6</sup>  
Carrier Material: Aluminum, anodized

+ growth  
- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time [min]	10.0	20.0	30.0	40.0	50.0	60.0	70.0	80.0	90.0	100.0	+
Result 1	+	+	+	+	+	+	+	-	-	+	neg
2	+	-	-	+	+	-	+	-	-	-	-
3	+	-	-	+	+	-	-	-	-	-	

- estimated D-value: 8.4 [min]
- Model Behavior: not OK

*Carrier Material,*


*Aluminium; Air Sampler*

### Roughness

- P-P: 1436 nm
- RMS: 290 nm
- $R_a$ : 236 nm

### Wettability

- high



# Carrier Material, Model Behavior

Test Organism: B. stearothermophilus  
Initial Population: > 1.0 x 10<sup>6</sup>  
Carrier Material: Aluminum, Air Sampler

+ growth  
- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Expousure [min]	10.0	20.0	30.0	40.0	50.0	60.0	70.0	80.0	90.0	100.0	+
Result 1	+	+	+	+	+	+	+	-	-	-	neg
2	+	+	+	+	+	-	-	-	-	-	-
3	+	+	+	+	+	-	-	-	-	-	

- estimated D-value: 7.9 [min]
- Model Behavior: OK




## *Carrier Material, Aluminium; Filling Line*

### Roughness

- P-P: 1478 nm
- RMS: 132 nm
- $R_a$ : 94 nm

### Wettability

- high



# Carrier Material, Model Behavior

Test Organism: B. stearothermophilus

Initial Population: > 1.0 x 10<sup>3</sup>

Carrier Material: Aluminum, anodized

+ growth

- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Expousure [min]	10.0	20.0	30.0	40.0	50.0	60.0	70.0	80.0	90.0	100.0	+
Result 1	+	+	+	+	+	+	+	+	+	+	neg
2	+	+	+	+	+	+	+	+	+	+	-
3	+	+	+	+	+	+	+	+	+	+	

- estimated D-value: > 33.7 [min]
- Model Behavior: not able to be evaluated

# *Selection of Bioindicators*

# *Carrier Material*

Carrier Material of BI	D-value Estimations [mins]	Model Behaviour
Glass	1.0 / 1.1	given
CrNi steel 1.4435, polished	1.3 / 0.9	given
CrNi steel 1.4301, unpolished	1.0 / 1.2	given
CrNi steel 1.4435, unpolished	1.0 / 1.4	given
CrNi steel 1.4301, polished	1.3 / 1.4	given
PVC, hard	1.0 / 1.8	given
PTFE	1.3 / 1.6	given
PE, UHMW	1.6 / 1.6	given
PP	1.3 / 2.0	given
PVC	2.0 / 1.6	given
Laminated foil 1, outside	1.6 / 2.5	given
PC	2.2 / 2.3	given
BI, commercially available	2.6 / 2.3	given
Tyvek	2.0 / 3.1	given
Laminated foil 2, outside	2.5 / 3.2	given
Butyl chaouchouc	2.9 / 3.1	given
Hypalon	3.0 / 4.1	given
HEPA-filter pad	3.6 / 3.6	given
PVC, soft	4.3 / 4.6	given
POM	4.6 / 4.4	given
Aluminium, anodized, commercially available	> 3.1 / 7.9	not given
Aluminium, anodized, air sampler	>8.3 / 10.1	given
Aluminium, anodized, filling line	>17.1 / > 33.7	not applicable

## *Selection of Bioindicators*

## *Carrier Material*



*Significantly higher* D-value for the various *aluminium* samples

## *Selection of Bioindicators*

## *Carrier Material*


### *Unsuitable Carrier Material*

- Negative influence on the *Resistance* of Bioindicators
- Negative influence on the *Model Behavior* of Bioindicators
- Materials are *spongy* or *porous*
- Test Organism *penetrate into* those Materials
- *not or only hard* to be reached by H<sub>2</sub>O<sub>2</sub> Decontamination


*Carrier Material      stainless steel, not polished*



*Carrier Material      stainless steel, not polished*




*Carrier Material      stainless steel, not polished*





*Carrier Material      stainless steel, not polished*



*Carrier Material*

*Aluminium; Air Sampler*




*Carrier Material*

*Aluminium; Air Sampler*



*Carrier Material*

*Aluminium; Air Sampler*




*Carrier Material*

*Aluminium; Air Sampler*




*Carrier Material*

*Aluminium anodized*




*Carrier Material*

*Aluminium anodized*



*Carrier Material*


*Aluminium anodized*






*Carrier Material*

*Aluminium anodized*



*Carrier Material*

*Aluminium anodized*



## *Selection of Bioindicators*

## *Carrier Material*

### *Unsuitable Carrier Material*


- Negative influence on the *Resistance* of Bioindicators
- Negative influence on the *Model Behavior* of Bioindicators
- Materials are *spongy* or *porous*
- Test Organism *penetrate into* those Materials
- *not or only hard* to be reached by H<sub>2</sub>O<sub>2</sub> Decontamination

The *Penetration* of the Decontamination Effect of the H<sub>2</sub>O<sub>2</sub> Process is

- *not ensured*
- *but also not expected*

# Selection of Bioindicators

# Carrier Material



# *Selection of Bioindicators Primary Packaging*

*Standards EN/ISO 14161*

„Primary Packaging should *not effect* the Inactivation and *protect* the inoculated Carrier against *Destruction* and *Contamination*“

## *H<sub>2</sub>O<sub>2</sub> Decontamination*

- semi permeable Membrane Tyvek
- permeable for H<sub>2</sub>O<sub>2</sub>
- Barrier for Contamination

# Primary Packaging

Test Organism: B. stearothermophilus  
Initial Population:  $> 1.0 \times 10^6$   
Carrier Material: Stainless Steel

+ growth  
- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time [min]	2.0	4.0	6.0	8.0	10.0	12.0	14.0	16.0	18.0	20.0	+
Result 1	+	+	+	+	-	-	-	-	-	-	neg
2	+	+	+	-	-	-	-	-	-	-	-
3	+	+	+	-	-	-	-	-	-	-	

- Primary Packaging: Tyvek
- estimated D-value: 1.2 [min]
- Model Behavior: OK

# Primary Packaging

Test Organism: B. stearothermophilus  
Initial Population: > 1.0 x 10<sup>6</sup>  
Carrier Material: Stainless Steel

+ growth  
- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time [min]	2.0	4.0	6.0	8.0	10.0	12.0	14.0	16.0	18.0	20.0	+
Result 1	+	+	-	-	-	-	-	-	-	-	neg
2	+	+	-	-	-	-	-	-	-	-	-
3	+	+	-	-	-	-	-	-	-	-	

- Primary Packaging: without primary Packaging
- estimated D-value: 0.8 [min]
- Model Behavior: OK

# *Selection of Biological Indicators*

*BI Composition has to reflect the Process Requirements*

- Test organisms      *Bacteria Spores* generally suitable
- Initial Population      *Preparation* has to be suitable for Process
- Carrier Material      *Properties* has to be suitable for Process
- Primary Packaging      *Properties* has to be suitable for Process
  
- Bioindicators      *reflect* the Decontamination effect *realistic*
- Bioindicators      show a *conform* Model Behavior



## *Selection of Biological Indicators*

- Bioindicators show *unrealistic* high Resistances
- Bioindicators show *Artifacts* in Model Behavior

## *Interpretation of the Process*

- Bioindicators show *Process- Performance und Boundaries*
- Bioindicators are *to high Process Challenge*
- Bioindicators show *Process to be overrated*
- Bioindicators show *Kind and Extent of acceptable Bioburden*

*Definition: Process- Performance und Expectations*

# *Biological Indicator for gaseous H<sub>2</sub>O<sub>2</sub>*

Testorganism: B. stearothermophilus ATCC 12980, min 4.5x10<sup>5</sup>  
Carrier / Package: Glasfibre / Tyvek  
Spezified D-value: 1.5 [min]

+ growth  
- no growth

## Run 01

exposure		01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]		5.0	10.0	15.0	20.0	25.0	30.0	35.0	40.0	45.0	50.0	+
result	1	+	+	+	+	-	+	+	-	-	-	neg
	2	+	+	-	-	-	-	-	-	-	-	-
	3	+	-	-	-	-	-	-	-	-	-	

## Run 02

exposure		01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]		10.0	20.0	30.0	40.0	50.0	60.0	70.0	80.0	90.0	100.0	+
result	1	+	-	+	-	+	-	-	-	-	-	neg
	2	+	-	-	-	-	-	-	-	-	-	-
	3	+	-	-	-	-	-	-	-	-	-	

# *Biological Indicator for gaseous H<sub>2</sub>O<sub>2</sub>*

Testorganism: B. stearothermophilus ATCC 7953, min 1.0 x10<sup>6</sup>  
 Carrier / Package: Glasfibre / Tyvek, PE  
 Spezified D-value: 3.1 [min]

+ growth  
 - no growth

## Run 01

exposure		01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]		6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
result	1	+	-	+	+	+	+	+	-	-	+	neg
	2	+	-	+	+	+	+	+	-	-	+	-
	3	-	-	-	+	-	+	-	-	-	+	

## Run 02

exposure		01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]		10.0	20.0	30.0	40.0	50.0	60.0	70.0	80.0	90.0	100.0	+
result	1	+	+	-	+	+	+	+	-			neg
	2	+	+	-	-	+	-	-	-			-
	3	+	-	-	-	-	-	-	-			

# *Selection of Biological Indicators      Steam*

Test Organism:      B. stearothermophilus

Initial Population:      > 1.0 x 10<sup>6</sup>

Carrier Material:      Paper

+ growth

- no growth

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time [min]	10.0	20.0	30.0	40.0	50.0	60.0	70.0	80.0	90.0	100.0	+
Result 1	+	+	+	+	+	+	+	+	+	+	neg
2	+	+	+	+	+	+	+	+	+	+	-
3	+	+	+	+	+	+	+	-	-	+	

- estimated D-value:      > 17.0 [min]
- Model Behavior:      not able to be evaluated

# *Selection of Biological Indicators*

Test Organism: B. stearothermophilus  
 Initial Population: > 1.0 x 10<sup>6</sup>  
 Carrier Material: CrNi- Stahl  
 Primary packaging: Tyvek

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time [min]	6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
Result 1	+	+	+	+	-	-	-	-	-	-	neg
2	+	+	+	-	-	-	-	-	-	-	-
3	+	+	-	-	-	-	-	-	-	-	

- estimated D-value: = 2.0 [min] + growth
- Model Behavior: = OK - no growth

# *Selection of Biological Indicators*

## *Alternative Sterilization Methods*

- Resistance and Model Behavior *describes* Bioindicator
- Process- Performance and Expectations *defines* Bioindicator
- careful Selection of *Composition* Bioindicators

*Artifacts* in the Results of Bioindicators *Process Boundaries*

# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination of Isolator Systems*



Volker Sigwarth  
Skan AG, Switzerland



Alexandra Stärk  
Novartis Pharma AG, Switzerland

# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination of Isolator Systems*

Volker Sigwarth  
Skan AG  
Switzerland




## *Process Development*



# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination Process Development*

- Process Development of alternative Sterilization Method



# *Investigation on Germ Reducing Parameters of the H<sub>2</sub>O<sub>2</sub> Decontamination Method*

- *Process Development*
- Establish parameters of microbial reducing effects of H<sub>2</sub>O<sub>2</sub> decontamination
- Target value is the *D-value* of a defined microbiological system
- Using the *Design of Experiment* method as statistical tool
- Statistical significant quantification of the influence of effects
- Correlation of *Process Parameter* versus *Kill Effect*


# *Investigation on Germ Reducing Parameters of the H<sub>2</sub>O<sub>2</sub> Decontamination Method*

- Comprehension of the decontamination process
- The *nice and need to have* parameters for the decontamination success
- Hints for Design, Qualification and Monitoring
- Method to *describe, develop and quantify* H<sub>2</sub>O<sub>2</sub> decontamination cycles

# *Investigation on Germ Reducing Parameters of the H<sub>2</sub>O<sub>2</sub> Decontamination Method*

- Test - System and - Equipment
- Microbiological System
- Decontamination Process
- Selection of the experimental Factors
- Design of Experiment
- Result of Investigation and Interpretation
- Summary

# Test - Isolator



Material: stainless steel, glass

Volume: 1,4 m<sup>3</sup> (40 ft<sup>3</sup>)

# *Sensors*

## *Process*

- Temperature [°C]
- Humidity [% rH]
- Air Velocity [m/s]
- Pressure [Pa]
- Mass / Balance [g]

## *H<sub>2</sub>O<sub>2</sub> Gas-Concentration*

- Electro-chemical Sensor
- UV-Spectrometer
- IMS-Spectrometer
- NIR-Spectrometer
- Wet-chemistry Method


# *Microbiological System*

- Commercially available Biological Indicators
- *Bacillus stearothermophilus*      ATCC 12980
- Stainless Steel Carriers
- Tyvek Pouches

# Microbiological System

## Model of Microbial Reduction

- Initial Population [log-scale]
- Inactivation Time [min]
- Survival Curve
- D-value [min]
- Survival - Kill Window [min]





## *Minimized LSKM, Reactive Pattern*

exposure	01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]	6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
result 1	+	+	+	+	-	-	-	-	-	-	neg
2	+	+	+	-	-	-	-	-	-	-	-
3	+	+	-	-	-	-	-	-	-	-	

estimated D-value = 2.0 [min]

+ growth  
- no growth

## Limited Spearman Karber Method LSKM


exposure		01	02	03	04	05	06	07	08	09	10	
exposure time [min]		6.0	8.5	11.0	13.5	16.0	18.5	21.0	23.5	26.0	28.5	
result	1	+	+	+	+	+	+	-	-	-	-	neg
	2	+	+	+	+	+	+	-	-	-	-	-
	3	+	+	+	+	-	-	-	-	-	-	-
	4	+	+	+	+	-	-	-	-	-	-	pos
	5	+	+	+	-	-	-	-	-	-	-	+
	6	+	+	+	-	-	-	-	-	-	-	+
	7	+	+	+	-	-	-	-	-	-	-	
	8	+	+	+	-	-	-	-	-	-	-	
	9	+	+	-	-	-	-	-	-	-	-	
	10	+	-	-	-	-	-	-	-	-	-	

D-value = 2.06 [min]


95% CI D-value = 2.06 ± 0.22 [min]

+ growth  
- no growth

# *Decontamination Cycle*




# *Decontamination Cycle*




Phase 1:            Pre-conditioning  
                         to establish the initial conditions in the chamber

# *Decontamination Cycle*




Phase 2:            Conditioning  
                          to establish the decontamination effect

# *Decontamination Cycle*




Phase 3:            **Decontamination**  
to maintain stable decontamination effect  
to ensure the total bacterial reduction over time

# *Decontamination Cycle*




Phase 4:            **Aeration**  
to reach minimal residual H<sub>2</sub>O<sub>2</sub> concentration

# *Selection of Factors*





# *Selection of Factors*




## **Pre-conditioning**

Initial conditions in the chamber

*Temperature* [°C]

*Humidity* [%rH]

# *Selection of Factors*




**Conditioning**

Vaporized initial quantity of *pure* H<sub>2</sub>O<sub>2</sub>

*Quantity*

*[g/m<sup>3</sup>]*

# *Selection of Factors*



Decontamination

Rate of continuously vaporized *pure* H<sub>2</sub>O<sub>2</sub>     *Redose*     [%/h]

## *Selection of Factors*

Special

Concentration of  $H_2O_2$  Solution

$H_2O_2$  [%]

5 factors selected for the investigation

*A: Quantity of pure  $H_2O_2$  [g/m<sup>3</sup>]*

*B: Rate of Redose [%A/h]*

*C: Temperature [°C]*

*D: Humidity [%rH]*

*E: Concentration of  $H_2O_2$  [%]*


## *Range of Factors*

No:	Description	Unit	SP-	FF-	CP	FF+	SP+
A:	Quantity of pure H <sub>2</sub> O <sub>2</sub>	g/m <sup>3</sup>	4	5	6.5	8	9
B:	Rate of Redose	% A/h	20	40	70	100	120
C:	Temperature	°C	26	30	35	40	44
D:	Humidity	% rH	6	10	15	20	24
E:	Concentration of H <sub>2</sub> O <sub>2</sub>	%	30	35	42.5	50	55

# *Design of Experiment DoE*

- Fractional Factorial Plan
- Centre Point
- Star Point

FF  
CP  
SP




Independent Estimation of the Parameters for:

- Main Effects
- Quadratic Effects
- Interactions
- + 4 independent determinations of Centre Point

## *Range of Factors*

No:	Description	Unit	SP-	FF-	CP	FF+	SP+
A:	Quantity of pure H <sub>2</sub> O <sub>2</sub>	g/m <sup>3</sup>	4	5	6.5	8	9
B:	Rate of Redose	% A/h	20	40	70	100	120
C:	Temperature	°C	26	30	35	40	44
D:	Humidity	% rH	6	10	15	20	24
E:	Concentration of H <sub>2</sub> O <sub>2</sub>	%	30	35	42.5	50	55


# *Test Handling*



- place BI`s in the isolator      gastight wrapped




# *Test Handling*




- establish required initial conditions
- start vaporizing H<sub>2</sub>O<sub>2</sub> up to required quantity

# *Test Handling*



- expose BI`s to the inactivation atmosphere
- start redosing

# Test Handling



- remove BI`s out of the chamber in constant time intervals

# *Results and Interpretation*      *Centre Point*


## *Independent Centre Point Determinations*

- full LSKM      10 groups, 10 BI`s per group

<i>Run No.:</i>	<i>D-value [min]</i>
1	<i>1.75</i>
2	<i>1.80</i>
3	<i>1.74</i>
4	<i>1.71</i>


- *Constant* parameters lead to a *reproducible* decontamination effect

# Results and Interpretation Model




# Results and Interpretation

# Quantity (A)





# *Results and Interpretation*      *Quantity (A)*

## Plot of D-value versus Quantity



# Results and Interpretation

# Quantity (A)





- Quantity of vaporized  $H_2O_2$  is the *Main Effect* in decontamination process
- Nonlinear effect                      quadratic fit within range of model
- Asymptotic character                observed with larger range than model




# Results and Interpretation

# Quantity (A)



- Starting with constant initial conditions  
boost the quantity of vaporized  $H_2O_2$   
leads to a *Steady State* in decontamination effect


# Results and Interpretation Redose (B)




- The *higher* the rates of Redose (B) the *better* the decontamination effect

# Results and Interpretation

# Interaction (AB)




# *Results and Interpretation*      *Interaction (AB)*



- Redose (B) loses its influence at higher values of factor Quantity (A)
- Only high rates of Redose (B) improve a low decontamination effect


# Results and Interpretation      Interaction (AB)



- *Stability* of decontamination effect over time depends on the *Rate of Redose*


# Results and Interpretation

Temperature (C)  
Humidity (D)



# Results and Interpretation


# Temperature (C)



- The *lower* the Temperature (C) the *better* the decontamination effect
- Minor effect: change in D-value of only + 0.5 [min] within a temperature range of 10 [°C]

## Results and Interpretation

## Temperature (C)




- *Worst Case* positions for the decontamination effect are positions with *higher temperature*






# Results and Interpretation      Humidity (D)




- *Worst Case* positions for the decontamination effect are positions with *lower humidity*

# Results and Interpretation $H_2O_2$ (E)




# Results and Interpretation $H_2O_2$ (E)



- The *lower* the Concentration of  $H_2O_2$  Solution (E) the *better* the decontamination effect
- Minor effect: change in D-value of only + 0.6 [min] within a concentration range of 15 [%]

# Results and Interpretation $H_2O_2$ (E)



- Vaporization of the same quantity of pure  $H_2O_2$  produces a *better* decontamination *effect* at *lower concentrations* of the  $H_2O_2$  solution

## *Results and Interpretation*

One main *Interrelationship* between *all Factors* ?

A: Quantity	Steady State	higher <i>saturation</i> of gaseous phase
B: Redose	Stability	higher <i>saturation</i> of gaseous phase
C: Temp.	lower temperature	higher <i>saturation</i> of gaseous phase
D: Humidity	higher humidity	higher <i>saturation</i> of gaseous phase
E: H <sub>2</sub> O <sub>2</sub>	lower concentration	higher <i>saturation</i> of gaseous phase


- Decontamination effect depends on *Saturation of gaseous phase*
- *“Physical Pressure” from gaseous phase to surface*

## *Results and Interpretation*      *H<sub>2</sub>O<sub>2</sub> Gas-Concentration*

- No useful correlation to the microbial reduction
- Depends highly on the temperature
- Most sensors show the same shape of measurement curve,  
but highly differences in values

# Results and Interpretation

# $H_2O_2$ Gas-Concentration





## *Results and Interpretation*      *H<sub>2</sub>O<sub>2</sub> Gas-Concentration*

- No useful correlation to the microbial reduction
- Depends highly on the temperature
- Most sensors show the same shape of measurement curve, but highly differences in values
- *Calibration method !!!*
- *Not yet a tool for describe microbial reduction*

# *Results and Interpretation*      *H<sub>2</sub>O<sub>2</sub> Gas-Concentration*

## Measurement of process concentration

- useful as indicative measurement

## Measurement of residual H<sub>2</sub>O<sub>2</sub> concentration

- important for the final release of the decontaminated area

## *Description of Decontamination Effect*

The decontamination effect is described by:

- the results of the minimized LSKM
- the only useful sensor is the BI
- the minimized LSKM leads to quantifiable results
  
- the suitability of the BI has to be tested prior
- to “calibrate” the BI as

*Sensor for the decontamination effect*

## *Development and Quantification of H<sub>2</sub>O<sub>2</sub> Decontamination Cycles*

# *Selection of Biological Indicators*


Test Organism: B. stearothermophilus  
 Initial Population: > 1.0 x 10<sup>6</sup>  
 Carrier Material: CrNi- Stahl  
 Primary packaging: Tyvek

Group	01	02	03	04	05	06	07	08	09	10	pos
Exposure Time [min]	6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
Result 1	+	+	+	+	-	-	-	-	-	-	neg
2	+	+	+	-	-	-	-	-	-	-	-
3	+	+	-	-	-	-	-	-	-	-	

- estimated D-value: = 2.0 [min] + growth
- Model Behavior: = OK - no growth


# *Description of Cycle*

# *Physical*



# *Description of Cycle*

# *Physical*



Reproducibility of Decontamination Effect


*Initial conditions in chamber*

*Temperature* ± 5 [°C]

*Humidity* ± 5 [%rH]

# *Description of Cycle*

# *Physical*




Achieved Decontamination Effect

*Quantity of initial vaporized H<sub>2</sub>O<sub>2</sub>*

mass control, balance [g]

# Description of Cycle

# Physical



Stability of Decontamination Effect


Rate of Redose H<sub>2</sub>O<sub>2</sub>

mass control, balance [g/t]



# *Description of Cycle*

# *Physical*



Residual Concentration of H<sub>2</sub>O<sub>2</sub>

*Aeration curve*


low concentration sensor [ppm/t]

## *Description of Cycle*      *Decontamination effect*

- Based on a well known and defined BI
- Using the minimized LSKM as a methodical tool
- Proof the decontamination effect over the whole cycle
- Define the required values for cycle parameters


# *Description of Cycle*

# *Decontamination effect*



## *Description of Cycle*

## *Decontamination effect*




### Reproducibility of Decontamination Effect

*Physically defined by:*

*Temperature*  $\pm 5$  [°C]


*Humidity*  $\pm 5$  [%rH]

# Temperature and Humidity Mapping



# Description of Cycle

# Decontamination effect



## Achieved Decontamination Effect

*achieved D-value*

*defined*

*initial Quantity of  $H_2O_2$  [g]*

# *Achieved Decontamination effect*

Achieved decontamination effect

Quantity 5 [g/m<sup>3</sup>]

+ growth

- no growth

exposure	01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]	6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
result 1	+	+	+	+	+	+	+	+	-	-	neg
2	+	+	+	+	+	+	-	-	-	-	-
3	+	+	+	+	-	-	-	-	-	-	

# *Achieved Decontamination effect*

Achieved decontamination effect

Quantity 7.5 [g/m<sup>3</sup>]

+ growth


- no growth

exposure	01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]	6.0	9.0	12.0	15.0	18.0	21.0	24.0	27.0	30.0	33.0	+
result 1	+	+	+	-	-	-	-	-	-	-	neg
2	+	+	-	-	-	-	-	-	-	-	-
3	+	-	-	-	-	-	-	-	-	-	




# *Achieved Decontamination effect*

## Plot of D-value versus Quantity



## *Description of Cycle*

## *Decontamination effect*



### Stability of Decontamination Effect

*Stability of D-values*

*defined*

*Rate of Redose [g/t]*

# Stability of Decontamination effect

Rate of Redose 25 [%A/h]

+ growth  
- no growth

LSKM 1, exposition 5 [min] after conditioning

exposure	01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]	7.5	10.0	12.5	15.0	17.5						+
result 1	+	-	-	-	-						neg
2	+	-	-	-	-						-
3	-	-	-	-	-						

LSKM 2, exposition 30 [min] after conditioning

exposure	01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]	7.5	10.0	12.5	15.0	17.5						+
result 1	+	+	+	+	+						neg
2	+	+	+	+	+						-
3	+	+	+	+	+						

# Stability of Decontamination effect

Rate of Redose 100 [%A/h]

+ growth  
- no growth

LSKM 1, exposition 5 [min] after conditioning

exposure	01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]	7.5	10.0	12.5	15.0	17.5						+
result 1	+	-	-	-	-						neg
2	+	-	-	-	-						-
3	+	-	-	-	-						

LSKM 2, exposition 30 [min] after conditioning

exposure	01	02	03	04	05	06	07	08	09	10	pos
exposure time [min]	7.5	10.0	12.5	15.0	17.5						+
result 1	+	-	-	-	-						neg
2	+	-	-	-	-						-
3	-	-	-	-	-						

# *Description of Cycle*      *Worst Case Positions*


Definition of Worst Case Positions, (BI- Positions)

- |                                      |               |
|--------------------------------------|---------------|
| • Air Flow and Distribution          | Smoke Tests   |
| • Temperature, Humidity Distribution | Deco-Process  |
| • Chemo Indicator Mapping            | qualitative   |
| • deeply consider the Process        | Risk Analysis |

*Extreme Point of Application*      *no Assumptions*

# *Description of Cycle*

# *Decontamination effect*



## Duration of Decontamination Phase

*Worst Case Study*

*defined*

*D-value Worst Case*

*Duration of Decontamination [t]*

## *Description of Cycle*


## *Worst Case Positions*

### Assessment of Worst Case Positions

- 3x BI`s per defined WC Position      fractional results possible
- Deco time 10 x D-value best place      shortest cycle to total kill
- Use resulting BI pattern to estimate kill      - - - 10log, - - + 6log, + + + ??
- Calculate D-value worst case      Deco time / achieved reduction
- longer Deco time if required      10 x D-value worst case

# *Description of Cycle*

# *Decontamination effect*



Required log Reduction

*observed D-values*

*defined*

*achieved log Reduction*



# *H<sub>2</sub>O<sub>2</sub> Decontamination Process*

- Reproducible and stable process *if well developed*
- Development contents two systems
  1. Decontamination System
  2. Microbiological System
- Recognize the suitability of the microbiological system before hand
- Develop the decontamination cycle using:
  - the “calibrated BI *and*
  - the whole physical range *no assumption*
- Transparent and systematic process development and comprehension
- Strong evidence for the following step    Performance Qualification

# *Development and Quantification of H<sub>2</sub>O<sub>2</sub> Decontamination Cycles*

- **Validation Article**                      **PDA Journal**  
“Development and Quantification of H<sub>2</sub>O<sub>2</sub> Decontamination Cycles”  
Volker Sigwarth, Dr. Claude Moirandat, Skan AG  
PDA Journal, Vol. 54, July / August 2000

## *Reference Studies for H<sub>2</sub>O<sub>2</sub> Decontamination*

- Suitability study for all commercially available BI`s
- D-value studies for a wide range of spores and vegetative germs
- Customized determination of virus decontamination
- Special studies for devices and equipment
  
- Effect of Carrier Materials on the Resistances of BI`s
- Co-operation with Novartis, Stein


# *Material Study*

- Materials used in isolator as construction materials  
as disposables in routine work
- Inoculate with min  $1.0 \times 10^6$  spores of *B. stearothermophilus*
- Determination of D- value and Reactive Pattern
- Compare the results with commercially available BI
- Transfer results to different Production Isolator Systems
  
- Conclusion “HOW to handle the Material Question”

# Materials Study

# Results

Carrier Material of BI	D-value Estimations [mins]	Model Behaviour
Glass	1.0 / 1.1	given
CrNi steel 1.4435, polished	1.3 / 0.9	given
CrNi steel 1.4301, unpolished	1.0 / 1.2	given
CrNi steel 1.4435, unpolished	1.0 / 1.4	given
CrNi steel 1.4301, polished	1.3 / 1.4	given
PVC, hard	1.0 / 1.8	given
PTFE	1.3 / 1.6	given
PE, UHMW	1.6 / 1.6	given
PP	1.3 / 2.0	given
PVC	2.0 / 1.6	given
Laminated foil 1, outside	1.6 / 2.5	given
PC	2.2 / 2.3	given
BI, commercially available	2.6 / 2.3	given
Tyvek	2.0 / 3.1	given
Laminated foil 2, outside	2.5 / 3.2	given
Butyl chaouchouc	2.9 / 3.1	given
Hypalon	3.0 / 4.1	given
HEPA-filter pad	3.6 / 3.6	given
PVC, soft	4.3 / 4.6	given
POM	4.6 / 4.4	given
Aluminium, anodized, commercially available	> 3.1 / 7.9	not given
Aluminium, anodized, air sampler	>8.3 / 10.1	given
Aluminium, anodized, filling line	>17.1 / > 33.7	not applicable



*Significantly higher D-value for the various *aluminium* samples*

## *Materials Study*

## *Interpretation Aluminium*

- Extremely *porous surface* structure
- Suspension was *absorbed into surface*
- H<sub>2</sub>O<sub>2</sub> decontamination is *not able to penetrate into a surface*

➡ *No* or only a *bad* inactivation effect

➡ *Not suitable for H<sub>2</sub>O<sub>2</sub> decontamination*

*Carrier Material*


*Aluminium anodized*






*Carrier Material*

*Aluminium anodized*



# Materials Study

# Results



## *Materials Study*


## *Transferability*

• 1 x Reference Isolator	1,4 m <sup>3</sup>
• 2 x Sterility Test Isolator	2,2 m <sup>3</sup>
• 5 x Filling Isolator	9,0 m <sup>3</sup>
• 2 x Material Pass Through	15,0 m <sup>3</sup>

- IQ / OQ finished
- Cycle development based on commercial BI`s finished
- Initial condition and cycle parameter defined

# Materials Study

# Results



## *Materials Study*

- 1 x Reference Isolator
- 2 x Sterility Test Isolator
- 5 x Filling Isolator
- 2 x Material Pass Through


## *Transferability*

1,4 m<sup>3</sup>  
2,2 m<sup>3</sup>  
9,0 m<sup>3</sup>  
15,0 m<sup>3</sup>

	D-value [mins] Reference Isolator	D-value [min] Serility Test Isolator	D-value [min] H <sub>2</sub> O <sub>2</sub> Material Pass Through	D-value [min] Filling Isolator
B. stearothermophilus on Glass	1.45 ± 0.07	1.41 ± 0.04	0.96 ± 0.05	1.18 ± 0.20
Biological Indicator Reference	2.28 ± 0.06	2.22 ± 0.06	1.62 ± 0.08	1.72 ± 0.07
B. stearothermophilus on POM	4.33 ± 0.20	4.00 ± 0.30	3.38 ± 0.16	3.41 ± 0.23

# Materials Study


# Implementation



## *Materials Study*

## *Summary*

- Reproducible *different resistances* on *different carrier materials*
- Some Materials *seems* to be *not suitable* for H<sub>2</sub>O<sub>2</sub> decontamination



*borders of the process*

- Resistances were *transferable* to *different systems*
- Cycle development and qualification based on *commercial BI`s*
- *Knowledge* of all inactivation factors

# *Effect of Carrier Materials on the Resistance of B.stearothermophilus to gaseous H<sub>2</sub>O<sub>2</sub>*

- Research Article PDA Journal

“Effect of Carrier Materials on the Resistance of Spores of Bacillus Stearothermophilus to gaseous Hydrogen Peroxide”

Volker Sigwarth, Skan AG

Alexandra Stärk, Novartis Pharma AG

PDA Journal, Vol. 57, No.1 January / February 2003

- As Reference Study for Isolator Validation



# *Effect of Carrier Materials on the Resistance of B.stearothermophilus to gaseous H<sub>2</sub>O<sub>2</sub>*

- Implementation in PIC/S Recommendation

“The carrier type e.g. plastic, paper, metal or other, of the biological indicator organism should be relevant to the materials being gassed or shown to be irrelevant”

“If studies have been carried to show that lethality on carrier type a is similar to materials c, d, e, etc. with a similar sporicidal process, this would mean that *in house studies need not be carried out.*”

“The data would need to be from a *reputable* source”

## *Materials Study*

## *FDA; Comment*

- FDA Comment, R. Friedman, ISPE Washington Conference 2003

Reproducibly “different resistances of *B. stearothermophilus* to a H<sub>2</sub>O<sub>2</sub> decontamination on different carrier materials”.

“Because of their surface structure and properties, certain materials seems to be not suitable for the H<sub>2</sub>O<sub>2</sub> decontamination ...”

*But... resistances were found to be essentially “transferable”, so development work may greatly reduce the need to address extensively during decontamination cycle validation.*

Development of decontamination cycle parameters *should incorporate knowledge of material-effects*. Commercial BI used as control for this study.

Authors conclude that no single factor on its own is normally responsible for material effect

## *Development of alternative Sterilization Methods*

- Book Chapter            “Process Development of alternative Sterilization Methods”

Title:                            “Contamination Control in Parenteral Processing”

Published by:                Marcel Dekker Inc., USA

Edited by:                    Kevin Williams, Eli Lilly & Co., USA

Author:                        Volker Sigwarth, Skan AG, Switzerland

published:                    Middle of 2004

# *Decontamination of Isolator Systems*

## *Risks in working with Biological Indicators*

- “total kill approach”                      one surviving BI stops production
- Bioindicators are hand made              reliable resistance of single BI
- use enough Bioindicators                  You will find a survivor

*How to handle surviving Bioindicators during qualification work ?*

# *Decontamination of Isolator Systems*

## *Risks in working with Biological Indicators*

- “total kill approach”
- Where come the requirement “total kill” from ?
- Steam Sterilization *versus* Isolator Decontamination

## *Regulatory Requirements and Possibilities*

# *Decontamination of Isolator Systems*


## *Regulatory Requirements; FDA Aseptic Guideline*

“Normally a four- to six-log reduction can be justified depending on the application. The specific spore titer used and the selection of BI placement site should be justified. For example, demonstration of a four-log reduction should be sufficient for controlled , very low bioburden material introduced into a transfer isolator including wrapped sterile supplies that are briefly exposed to the surrounding environment”


# Microbiological System

## Model of Microbial Reduction

- Initial Population [log-scale]
- Inactivation Time [min]
- Survival Curve
- D-value [min]
- Survival - Kill Window [min]



# Microbiological System





# *Decontamination of Isolator Systems*

*Regulatory Requirements; USP 28, <1208>*

The sterilization methods used to treat isolators, test articles, and sterility testing supplies are capable of reproducibly yielding a six-log kill against an appropriate, highly resistant biological indicator (BI; see *Biological Indicators for Sterilization* (1035)), as verified by the fraction negative or total kill analysis methods. Total kill analysis studies are suitable for BIs with a population of  $10^4$  spores per unit, while fraction negative studies are suitable for BIs with a population of  $10^5$  or greater. A sufficient number of BIs are used to prove statistical reproducibility and adequate distribution of the sterilizing agent. Particular attention is given to areas that pose problems relative to the concentration of the agent. A larger number of BIs are used in isolators that are heavily loaded with equipment and materials. Also, when it is not possible to use one or more calibrated sensors to directly measure the concentration of the sterilizing agent, the placement of additional BIs is considered. The ability of the process to reproducibly deliver a six-log kill is confirmed in three consecutive validation studies.

# *Decontamination of Isolator Systems*

## *Possibilities in working with Biological Indicators*

- “fractional kill approach”
- relation between positive and negative Bioindicator results
- proof the required log reduction

## *Successfully applied for Filling and Sterility Test Isolators*

# *Decontamination of Isolator Systems*

## *Problems in using “fractional kill approach”*

- Basics of microbial inactivation has to be deeply understood
- Argumentations during audits and with authorities often much harder
- Steam Sterilization *versus* Isolator Decontamination
- Root cause of positive result ?    BI, Cycle Parameter, Isolator System

*Will be the future approach for isolator cycle qualification*

# *Decontamination of Isolator Systems*

## *Possibilities in working with Biological Indicators*

- Initial control of Bioindicators                      BI Quality Control
- D-value based Cycle Development                      quantifiable for each phase
- one unexpected positive BI result at one position
- followed by one more run with multiple BI samples at this position

*Successfully applied "Back Up Tool" for Cycle Qualification*

# *Biological Indicators and H<sub>2</sub>O<sub>2</sub> Decontamination of Isolator Systems*

Volker Sigwarth  
Skan AG  
Switzerland



## *Process Development*